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Silver

Qualitative analysis

10 ml of the water sample is treated in a test-tube with a few drops of Feigl's reagent, 5 ml diethyl ether superimposed, and the mixture shaken. In the presence of silver ions a reddish-brown film forms at the interface between the two phases. Lower limit of detection: about 0.4mg/l.

Reagent preparation

Feigl's reagent: Dissolve 0.03g 5-(4-Dimethylaminobenzylidene) rhodanine GR in 100 ml Acetone GR.

Reagents

Cat. No.	Packages
14 Acetone GR	1l, 2 $\frac{1}{2}$ l
921 Diethyl ether GR	1l, 5l
3059 5-(4-Dimethylaminobenzylidene) rhodanine GR (reagent for silver).....	5g

Quantitative analysis

Photometrically with dithizone

Suitable for determination of silver ions in concentrations from 0.05 to 1.50mg Ag^+/l .

Interferences and pre-treatment

In order to achieve satisfactory results of measurement, all the apparatus should be cleaned extremely carefully before the determination by subjecting it to the action of Extran solution for one to two hours, and then rinsing successively with deionised water, dithizone solution, carbon tetrachloride and double-distilled water. To test for cleanness, 20 ml double-distilled water is introduced into the separating funnels before each determination, adjusted to pH about 9.5 with 1-2 drops of ammonia solution (0.91), and 0.5 ml dithizone solution added. After vigorous shaking for one minute the carbon tetrachloride phase should be colourless.

Interferences by zinc, cadmium and nickel ions in higher concentrations ($>100 \text{ mg/l}$) and by copper ions are eliminated by addition of 2 ml Titriplex III solution per sample volume.

In the presence of sulfide ions the sample is evaporated to dryness on a water bath with 1 ml nitric acid (about 1.40). The residue is moistened with 2 ml nitric acid (about

1.40), evaporated to dryness once again, and taken up in 50 ml doubledistilled water, with heating.

The interference of organic substances and iron (III) ions is excluded by addition of hydroxylammonium chloride solution; this method is not applicable to solutions acidified with nitric acid. Larger amounts of organic constituents (potassium permanganate-consuming capacity >25 mg/l) are decomposed by evaporation with sulfuric acid (about 1.55) and a few drops of nitric acid (about 1.40).

Procedure

Extreme care, and cleanness of apparatus, is necessary to avoid interferences; the determinations should always be carried out with the same cleaned apparatus.

Treat 100 ml test water in a 250 ml separating funnel with 15 ml nitric acid (about 1.035) or sulfuric acid (about 1.05), then with 50.0 ml dithizone solution, and shake vigorously for 5 min. If the carbon tetrachloride phase develops a golden yellow colour, a smaller starting volume of the test water should be taken. After complete separation of the phases, run the organic extract into a second separating funnel, treat with 10 ml ammonia solution (about 0.98), and shake for 2 min in a shaking machine. Dry the outlet tube of the separating funnel with a small roll of filter paper. Then allow some of the organic phase to run out in order to rinse the outlet tube. Filter the carbon tetrachloride layer through a plaited filter into a 1-cm cell, and measure at 470 nm (e.g., in an ELKO II. or III with Filter S 47 E) against a blank treated in the same manner.

The reference curve is prepared from standard solutions of silver nitrate treated in the same manner as the water sample.

Calculation

The starting volume of the test water is to be taken into account in calculating the silver ion content in mg Ag^+ /l.

1 meq silver ions $\triangleq 107.88$ mg Ag^+

Reporting of results

The values are rounded off to 0.01 mg/l for a silver ion content of less than 1 mg/l, and to 0.1 mg/l for 1 mg/l or more.

Reagent preparations

All the solutions are prepared with double-distilled water.

Ammonia solution (about 0.98): Mix 110 ml Ammonia solution min. 25% (0.91) Suprapur with 400 ml double-distilled water.

Dithizone solution: Dissolve about 12–13 mg Dithizone (1,5-diphenylthiocarbazon) GR in 100 ml Carbon tetrachloride GR (for determinations with dithizone). —Dilute 10 ml of this stock solution with 100 ml Carbon tetrachloride GR (for determinations with dithizone). Only stable for about 1 day! —Store the solutions in brown bottles.

Extran solution: 2% aqueous solution of Extran MA 01 alkaline.

Hydroxylammonium chloride solution: Dissolve 20g Hydroxylammonium chloride GR in 100 ml double-distilled water. The solution is stable for about 1 week. Nitric acid (about 1.035): Mix 580 ml double-distilled water with 50 ml Nitric acid min. 65% (about 1.40) tested by the dithizone method GR.

Sulfuric acid (about 1.55): Carefully pour 350 ml Sulfuric acid 95-97% (about 1.84) GR into 310 ml double-distilled water.

Sulfuric acid (about 1.05): Mix 440 ml double-distilled water with 22 ml Sulfuric acid 95-97% (about 1.84) GR.

Silver standard solution: Treat 1.5746g Silver nitrate GR with 10 ml nitric acid (1.035), and make up to 1000 ml with double-distilled water. 1 ml \triangleq 1 mg Ag⁺. or:

Treat 1.4453 g Silver sulfate GR with 10 ml sulfuric acid (1.05), and make up to 1000 ml with double-distilled water. 1 ml \triangleq 1mg Ag⁺.

The dilution solutions necessary for plotting of the reference curve are freshly prepared from one of these solutions.

Titriplex III solution: Dissolve 5 g Titriplex III GR in double-distilled water to make 100 ml.

Reagents

Cat. No.	Packages
9532 Alkalit [®] indicator strips pH 7.5-14.0 non-bleeding graduated in 0.5 pH units	1 pack
5428 Ammonia solution min. 25% (0.91) Suprapur [®]	250 ml, 1 l
3092 Dithizone (1,5-diphenylthiocarbazone GR	5 g, 25 g, 100 g
7555 Extran [®] MA 01 alkaline	1 l, 5 l
4616 Hydroxylammonium chloride GR	100 g, 250 g, 1 kg
456 Nitric acid min. 65% (about 1.40) tested by the dithizone method GR	1 l, 2 1/2 l
731 Sulfuric acid 95-97% (about 1.84) GR	1 l, 2 1/2 l
1512 Silver nitrate GR	25 g, 100 g, 250 g, 1 kg
1509 Silver sulfate GR	25 g, 100 g
2208 Carbon tetrachloride GR (for determinations with dithizone)	1 l, 2 1/2 l
8418 Titriplex [®] III [(ethylenedinitrilo) tetraacetic acid disodium salt] (dihydrate) GR	100 g, 250 g, 1 kg, 5 kg

Mercury

Quantitative analysis

Photometrically with dithizone

Suitable for determination of mercury ions in concentrations between 0.01 and 2.0 mg Hg^{2+} /l.

Interferences and pre-treatment

Silver, gold, palladium and platinum ions interfere at all concentrations.

The interference of larger concentrations of copper ($\text{Cu} : \text{Hg} > 10^3$) is eliminated as described under Procedure.

For determination of the total mercury organically bound mercury must be mineralised before the quantitative estimation, and the volatility of most mercury compounds must be borne in mind: Boil 100 ml of the water sample for 2h under reflux with 10 ml sulfuric acid (about 1.84) and sufficient potassium permanganate solution to leave a slight violet colouration at the end of the oxidation process. After cooling, decolourise, or dissolve precipitated manganese oxide, with a slight excess of hydroxylammonium chloride solution, boil, and cool to room-temperature.

Since mercury ions have a great tendency to adsorption, all glass apparatus should be carefully rinsed with dilute nitric acid, and then deionised water, before use.

Procedure

Introduce 100 ml, or a smaller volume made up to 100 ml with deionised water, of the original water sample acidified with 1 N sulfuric acid to pH 0-1, or 100 ml of the sample pre-treated as necessary, into a separating funnel, and extract with two 10.0-ml portions of dithizone solution. The last extract must remain green; if not, a smaller volume of the water made up to 100 ml with deionised water should be used.

If the water sample contains larger amounts of copper ($\text{Cu} : \text{Hg} > 10^3$), extract fractionally in the separating funnel as often as is required for the pure violet colour of copper dithizonate to appear. Shake the combined organic extracts twice for 2 min in each case with 15 ml 1 N sulfuric acid and 1 ml potassium permanganate solution. The mercury goes into the aqueous phase, which is decolourised with a few drops of hydroxylammonium chloride solution. Now extract the aqueous phase with two 10.0 ml portions of dithizone solution.

Wash the combined organic extracts (from the treatment described in paras. 1 and/or 2) free from dithizone with three 10-ml portions of dilute ammonia solution, and then desensitise by shaking vigorously with 10 ml dilute acetic acid. After separation of the phases, drain the organic phase through a filter paper into a 2-cm cell, cover the cell immediately, and measure at 485 nm (e.g., in an ELKO II or III with Filter S 49 E)

against a blank treated in the same manner. The reference curve is prepared from standard solutions of mercury (II) chloride or dilution solutions of Mercury standard solution Titrisol treated in the same manner as the water sample.

Calculation

Any dilution of the water sample to 100 ml starting volume is to be taken into account in the calculation.

Reporting of results

The values are rounded off to 0.01 mg/l.

1 meq mercury(II) ions \triangleq 100.30 mg Hg²⁺.

Reagent preparations

Ammonia solution, dilute: Dilute 5 ml Ammonia solution min. 25% (0.91) GR to 1000 ml with deionised water.

Dithizone solution: Dissolve 0.02g Dithizone GR by shaking for a longish period in 100 ml Carbon tetrachloride GR (for determinations with dithizone), and then extract with four 50-ml portions of dilute ammonia solution.

Filter the combined aqueous extracts through a plaited filter, add 100 ml Carbon tetrachloride GR and 1 ml Hydrochloric acid min. 25% (about 1.125), and shake immediately. Wash the green organic phase with three 50 ml portions of deionised water, and filter again. The solution is stable for about 4 weeks.

For the analysis a dilute solution of dithizone is employed which should be freshly prepared before use in each case: Make 10 ml of the dithizone stock solution up to 200 ml with Carbon tetrachloride GR.

Acetic acid, dilute: Dilute 15 ml Acetic acid glacial min. 96% (about 1.06) GR to 100 ml with deionised water.

Hydroxylammonium chloride solution: Dissolve 20 g Hydroxylammonium chloride GR in deionised water to make 100 ml.

Potassium permanganate solution: Saturated solution at room temperature (about 6 g Potassium permanganate GR in 100 ml deionised water).

Mercury standard solution: Dissolve 1.354 g Mercury (II) chloride GR in 1 N Sulfuric acid to make 1000 ml. 1 ml \triangleq 1 mg Hg²⁺.

For preparation of the reference curve fresh dilutions of this solution are made with 1 N Sulfuric acid as required.

The Mercury standard solution Titrisol can also be used in a similar manner. Nitric acid, dilute: Mix 100 ml Nitric acid min. 65% (about 1.40) tested by the dithizone method GR with 300 ml deionised water.

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